

Anti-CLIC1 Rabbit antibody



Catalog Number: 167603

Product name

Rabbit IgG

Anti-CLIC1 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

CLIC1 Rabbit monoclonal antibody

ICC/IF: 1/20

Preparation

FC: 1/50

Antigen: Recombinant protein of human CLIC1

IP: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage


Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of CLIC1 in HeLa,CHO-K1,C6,Ramos using CLIC1 Rabbit mAb(1:1000 diluted)