

# GOSR2/Membrin antibody



Catalog Number: 107303

## Product name

GOSR2/Membrin antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Immunogen

[Human GOSR2 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

GOSR2/Membrin Mouse Monoclonal antibody.  
Positive WB detected in HeLa cells. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of GOSR2/Membrin recombinant protein (Accession Number: NM\_004287). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107303(GOSR2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:107303(GOSR2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human liver slide using Catalog

# GOSR2/Membrin antibody



Catalog Number: 107303

---

No:107303(GOSR2 Antibody) at dilution of 1:50