

PLSCR3 antibody



Catalog Number: 114013

Product name

ELISA, IHC

PLSCR3 antibody

Dilutions

Specificity

Recommended Dilution:


Human; other species not tested.

IHC: 1:20-1:200

Antibody description

PLSCR3 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human skeletal muscle tissue.

Validations




Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_020360). Purification method: Protein A purified.

Immunohistochemical of paraffin-embedded human kidney using Catalog No:114013(PLSCR3 antibody) at dilution of 1:50 (under 10x lens)

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human kidney using Catalog No:114013(PLSCR3 antibody) at dilution of 1:50 (under 40x lens)

Applications