

TFIIIE beta antibody



Catalog Number: 116013

Product name

TFIIIE beta antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TFIIIE beta Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A375 cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of TFIIIE beta recombinant protein (Accession Number: XM_024447138). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116013(TFIIIE beta antibody) at dilution of 1:800