

Anti-TGF alpha Rabbit antibody



Catalog Number: 167013

Product name

Anti-TGF alpha Rabbit antibody

WB: 1/1000

IHC: 1/500

Specificity

Human

ICC/IF: 1/20

IP: 1/20

Antibody description

TGF alpha Rabbit polyclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human TGF alpha

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP

Dilutions



Immunohistochemistry of TGF alpha in paraffin-embedded Human colon cancer tissue using TGF



alpha Rabbit pAb at dilution 1/50

Western blot detection of TGF alpha in K562 cell lysates using TGF alpha Rabbit pAb(1:1000 diluted). Predicted band size: 17kDa. Observed band size: 17kDa.