

# BAIAP2L1 antibody



Catalog Number: 108413

## Product name

BAIAP2L1 antibody

## Immunogen

[Human BAIAP2L1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

BAIAP2L1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of BAIAP2L1 recombinant protein (Accession Number: NM\_018842). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108413(BAIAP2L1 Antibody) at dilution of 1:600