

Anti-DAPK3 Mouse antibody



Catalog Number: 169113

Product name

Anti-DAPK3 Mouse antibody

Specificity

Human

Antibody description

DAPK3 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human DAPK3 (AA: 28-161) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

IHC-P, FC


Validations



Flow cytometric analysis of A431 cells using DAPK3 mouse mAb (green) and negative control



Western blot analysis using DAPK3 mAb against




human DAPK3 (AA)

Western blot analysis using DAPK3 mAb against HEK293 (1) and DAPK3 (AA)

Anti-DAPK3 Mouse antibody



Catalog Number: 169113



Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using DAPK3 mouse mAb with DAB staining.