

SOX17 antibody



Catalog Number: 115513

Product name

SOX17 antibody

Immunogen

[Human SOX17 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SOX17 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue, C6 cells, PC-3 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of SOX17 recombinant protein (Accession Number: BC111365). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:115513(SOX17 Antibody) at dilution of 1:600