

# Anti-Eph receptor B1/NET Rabbit antibody



Catalog Number: 167913

## Product name

Anti-Eph receptor B1/NET Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

Eph receptor B1/NET Rabbit monoclonal antibody

## Dilutions


WB: 1/1000

IP: 1/20

## Preparation

Antigen: A synthetic peptide of human Eph receptor B1/NET

## Validations



## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

Western blot detection of Eph receptor B1/NET in A549,MCF-7,U251 cell lysates using Eph receptor B1/NET Rabbit mAb(1:1000 diluted).Predicted band size:108kDa.Observed band size:135kDa.