

ALKBH8 Polyclonal Antibody



Catalog Number: 165213

Product name

ALKBH8 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ALKBH8

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ALKBH8 (NP_620130.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:100

Validations




Western blot - ALKBH8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ALKBH8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - ALKBH8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using ALKBH8 antibody at dilution of



1:100 (40x lens).


ALKBH8 Polyclonal Antibody



Catalog Number: 165213

Immunohistochemistry - ALKBH8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using ALKBH8 antibody at



dilution of 1:100 (40x lens).

Immunofluorescence - ALKBH8 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ALKBH8 antibody . Blue: DAPI for nuclear staining.