

PRPF31 Polyclonal Antibody



Catalog Number: 164613

Product name

PRPF31 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PRPF31

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PRPF31 (NP_056444.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - PRPF31 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRPF31 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - PRPF31 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using PRPF31 antibody at dilution of 1:100



(40x lens).


PRPF31 Polyclonal Antibody



Catalog Number: 164613

Immunohistochemistry - PRPF31 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using PRPF31 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - PRPF31 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using PRPF31 antibody at dilution of 1:100 (40x lens).