

UNC45A antibody



Catalog Number: 116613

Product name

UNC45A antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

UNC45A Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IHC detected in human cervix tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of UNC45A recombinant protein (Accession Number: XM_011521779). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116613(UNC45A antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human cervix using Catalog No:116613(UNC45A antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded

UNC45A antibody




Catalog Number: 116613

human cervix using Catalog No:116613(UNC45A antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using UNC45A antibody Catalog No:116613 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-UNC45A (IP:Catalog No:116613, 4ug; Detection:Catalog No:116613 1:1000) with HeLa cells lysate 1400ug.