

Anti-PTP4A2/PRL2 Rabbit antibody



Catalog Number: 167613

Product name

Anti-PTP4A2/PRL2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

PTP4A2/PRL2 Rabbit monoclonal antibody

Dilutions


WB: 1/1000

IP: 1/20

Preparation

Antigen: A synthetic peptide of human PTP4A2/PRL2

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of PTP4A2/PRL2 in HeLa,A549,HL-60 cell lysates using PTP4A2/PRL2 Rabbit mAb(1:1000 diluted).Predicted band size:19kDa.Observed band size:19kDa.