

LRAT antibody



Catalog Number: 112313

Product name

LRAT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LRAT Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, mouse testis tissue. Observed molecular weight by Western-blot: 26 kDa, 50-54 kDa

Preparation

This antibody was obtained by immunization of LRAT recombinant protein (Accession Number: NM_004744). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:112313(LRAT antibody) at dilution of 1:500