

Anti-COX5B Rabbit antibody



Catalog Number: 167313

Product name

Anti-COX5B Rabbit antibody

IHC: 1/200

ICC/IF: 1/200

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

Antibody description

COX5B Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human COX5B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC, IP


Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of COX5B in K562, Rat Brain, C6, 3T3, Hela cell lysates using COX5B Rabbit pAb (1:1000 diluted). Predicted band size: 14 kDa. Observed band size: 14 kDa.



Immunohistochemistry of COX5B in paraffin-embedded Human colon cancer tissue using COX5B Rabbit pAb at dilution 1/50