

Drosha antibody



Catalog Number: 110023

Product name

Polyclonal

Drosha antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

Drosha Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 90 kDa, 151 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_013235). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110023(Drosha antibody) at dilution of 1:800