

ZC3H12A antibody



Catalog Number: 116923

Product name

ZC3H12A antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

ZC3H12A Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in Raji cells, mouse spleen tissue. Positive IHC detected in human small intestine tissue. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of ZC3H12A recombinant protein (Accession Number: NM_025079). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116923(ZC3H12A Antibody) at dilution of 1:300



IP Result of anti-ZC3H12A (IP:Catalog No:116923, 4ug; Detection:Catalog No:116923 1:300) with Raji cells lysate 2000ug.



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog

ZC3H12A antibody



Catalog Number: 116923

No:116923(ZC3H12A Antibody) at dilution of 1:50 (under 40x lens)