

ENPP1 antibody



Catalog Number: 110323

Product name

ENPP1 antibody

Immunogen

[Human ENPP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

ENPP1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of ENPP1 recombinant protein (Accession Number: BC059375). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:110323(ENPP1 antibody) at dilution of 1:1000