

Product name

humanin antibody

Specificity

Human; other species not tested.

Antibody description

humanin Rabbit Polyclonal antibody. Positive IHC detected in human small intestine tissue, human colon tissue, human heart tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NC_012920).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

Validations


Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:111723(humanin Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:111723(humanin Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:111723(humanin Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:111723(humanin Antibody) at dilution of 1:50 (under 40x lens)