

STK38 antibody



Catalog Number: 115723

Product name

STK38 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

STK38 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HeLa cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007271).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115723(STK38 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HEK-293 cells using Catalog No:115723(STK38 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)