

Anti-phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit antibody



Catalog Number: 167723

Product name

Rabbit IgG

Anti-phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit antibody

Applications

WB, IHC-P

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/1000

Antibody description


phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit monoclonal antibody

IHC: 1/50

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser181 of human TAOK1

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of TAOK1 (phospho-Ser181)/TAOK2 (phospho-Ser181)/TAOK3 (phospho-Ser177) in HeLa,A549,HL-60,U2OS,C6 cell lysates using TAOK1 (phospho-Ser181)/TAOK2 (phospho-Ser181)/TAOK3 (phospho-Ser177) Rabbit mAb(1:1000 diluted).Predicted band size:105,116,138kDa.Observed band size:105,116,138kDa.