

Anti-Cytokeratin 19 Rabbit antibody



Catalog Number: 167433

Product name

Anti-Cytokeratin 19 Rabbit antibody

Dilutions

WB: 1/5000

Specificity

IHC: 1/20

Human, Mouse

ICC/IF: 1/20

Antibody description

Cytokeratin 19 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Cytokeratin 19

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF

Validations



Immunohistochemistry of Cytokeratin 19 in paraffin-embedded Human colon cancer tissue using Cytokeratin 19 Rabbit pAb at dilution 1/20



Western blot detection of Cytokeratin 19 in K562 cell lysates using Cytokeratin 19 Rabbit pAb(1:1000 diluted). Predicted band size: 44kDa. Observed band size: 41kDa.