

Dynamitin antibody



Catalog Number: 110133

Product name

Dynamitin antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dynamitin Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, MCF7 cells, mouse brain tissue, mouse lung tissue. Positive IHC detected in human gliomas tissue. Positive FC detected in HeLa cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of Dynamitin recombinant protein (Accession Number: NM_001261413). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110133(dynactin-2 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110133(dynactin-2 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human gliomas using Catalog

Dynamitin antibody



Catalog Number: 110133

No:110133(dynactin-2 antibody) at dilution of 1:50
(under 40x lens)



1X10⁶ HeLa cells were stained with 1ug dynactin-2 antibody (Catalog No:110133, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.