

KRT40 antibody



Catalog Number: 112133

Product name

KRT40 antibody

Immunogen

[Human KRT40 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

KRT40 Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of KRT40 recombinant protein (Accession Number: NM_001389244). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:112133(KRT40 Antibody) at dilution of 1:300