

CAP2 antibody



Catalog Number: 108933

Product name

CAP2 antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:500-1:5000

IHC: 1:20-1:200

Antibody description

CAP2 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human skeletal muscle tissue. Positive WB detected in human kidney tissue, HepG2 cells, mouse heart tissue, mouse kidney tissue, mouse skeletal muscle tissue, mouse testis tissue, PC-3 cells, rat testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 53 kDa

Validations



Preparation

This antibody was obtained by immunization of CAP2 recombinant protein (Accession Number: NM_006366). Purification method: Antigen affinity purified.

human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108933(CAP2 antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP


Dilutions




IP Result of anti-CAP2 (IP:Catalog No:108933, 4ug; Detection:Catalog No:108933 1:1000) with mouse testis tissue lysate 4000ug.

CAP2 antibody

Catalog Number: 108933



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108933(CAP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108933(CAP2 Antibody) at dilution of 1:50 (under 40x lens)