

# DHX38 antibody



Catalog Number: 109933

## Product name

DHX38 antibody

## Immunogen

[Human DHX38 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DHX38 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 141 kDa

## Preparation

This antibody was obtained by immunization of DHX38 recombinant protein (Accession Number: XM\_011523484). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109933(DHX38 antibody) at dilution of 1:1000



IP Result of anti-DHX38 (IP:Catalog No:109933, 4ug; Detection:Catalog No:109933 1:700) with HeLa cells lysate 1400ug.