

# Anti-TTF2 Rabbit antibody

Catalog Number: 167233

## Product name

Anti-TTF2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TTF2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human TTF2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000

IHC: 1/50


ICC/IF: 1/20

FC: 1/20

## Validations



Immunocytochemistry of TTF2(green) in HeLa cells using TTF2 Rabbit pAb at dilution 1/50, and



Western blot detection of TTF2 in K562,Rat Brain,C6 cell lysates using TTF2 Rabbit pAb(1:1000 diluted).Predicted band size:130kDa.Observed band size:130kDa.