

transgelin/SM22 antibody



Catalog Number: 116333

Product name

transgelin/SM22 antibody

Immunogen

[Human TAGLN Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

transgelin/SM22 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, human colon tissue. Positive IHC detected in human heart tissue, human colon cancer tissue. Observed molecular weight by Western-blot: 19-22kd

Preparation

This antibody was obtained by immunization of transgelin/SM22 recombinant protein (Accession Number: NM_003186). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116333(TAGLN Antibody) at dilution of 1:1000




Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116333(TAGLN Antibody) at dilution of 1:50 (under 10x lens)

transgelin/SM22 antibody



Catalog Number: 116333



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116333(TAGLN Antibody) at dilution of 1:50 (under 40x lens)