

POLR2H antibody



Catalog Number: 114043

Product name

POLR2H antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POLR2H Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in HL-60 cells, HeLa cells, MCF7 cells, mouse kidney tissue, mouse liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 17kd

Preparation

This antibody was obtained by immunization of POLR2H recombinant protein (Accession Number: NM_001389583). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114043(POLR2H antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114043(POLR2H antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114043(POLR2H antibody) at dilution of 1:100

POLR2H antibody



Catalog Number: 114043

(under 40x lens)