

MRPS16 antibody



Catalog Number: 112843

Product name

MRPS16 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRPS16 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HEK-293 cells, Jurkat cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 15 kDa

Preparation

This antibody was obtained by immunization of MRPS16 recombinant protein (Accession Number: NM_016065). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112843(MRPS16 antibody) at dilution of 1:400



IP Result of anti-MRPS16 (IP:Catalog No:112843, 4ug; Detection:Catalog No:112843 1:500) with HEK-293 cells lysate 1480ug.