

Human IL-15 / IL15 / Interleukin 15 (His Tag) recombinant protein



Catalog Number: 503143

General Information

Protein Construction

A DNA sequence encoding the human IL15 (NP_000576.1) (Asn18-Ser162) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Activity

Measured in a cell proliferation assay using MO7e human megakaryocytic leukemic cells.
The ED₅₀ for this effect is typically 10-40 ng/mL.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human IL15 consists of 130 amino acids and predicts a molecular mass of 14.9 KDa. It migrates as an approximately 15 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL-15 / IL15 / Interleukin 15 Protein (His Tag) SDS-PAGE