

FBXL15 antibody



Catalog Number: 110543

Product name

FBXL15 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

FBXL15 Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 33 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of FBXL15 recombinant protein (Accession Number: BC002912). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:110543(FBXL15 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality