

Anti-PEAK3 antibody



Catalog Number: 176643

Product name

Anti-PEAK3 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to PEAK3

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, ICC, IHC-P

Dilutions

WB: 1:500-1:1,000

ICC: 1:50-1:200

IHC-P: 1:50-1:200

Validations




Fig1: Western blot analysis of C19orf35 on recombinant protein using anti-C19orf35 antibody at 1/1,000 dilution.




Fig2: ICC staining C19orf35 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig3: ICC staining C19orf35 in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig4: ICC staining C19orf35 in PC-3M cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.





Fig5: Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-C19orf35 antibody. Counter stained with



hematoxylin.

Fig6: Immunohistochemical analysis of paraffin-embedded human lung tissue using anti-C19orf35 antibody. Counter stained with hematoxylin.





Fig7: Immunohistochemical analysis of paraffin-embedded human uterus tissue using anti-C19orf35 antibody. Counter stained with



hematoxylin.

Fig8: Immunohistochemical analysis of paraffin-embedded human breast tissue using anti-C19orf35 antibody. Counter stained with hematoxylin.