

NUMBL antibody



Catalog Number: 107343

Product name

NUMBL antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

NUMBL Mouse Monoclonal antibody. Positive IF detected in L02 cells. Positive IHC detected in human gliomas tissue. Positive WB detected in L02 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of NUMBL recombinant protein (Accession Number: BC001794). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107343(NUMBLIKE Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:107343(NUMBLIKE Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded

NUMBL antibody



Catalog Number: 107343

human gliomas tissue slide using Catalog No:107343(NUMBLIKE Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of L02 cells using Catalog No:107343(NUMBLIKE Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)