

# VGAT antibody



Catalog Number: 116743

## Product name

VGAT antibody

## Immunogen

[Human SLC32A1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VGAT Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of VGAT recombinant protein (Accession Number: NM\_080552). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116743(SLC32A1 antibody) at dilution of



1:1000

IP Result of anti-SLC32A1 (IP:Catalog No:116743, 4ug; Detection:Catalog No:116743 1:500) with mouse brain tissue lysate 4000ug.