

# SAMD4B antibody



Catalog Number: 115053

## Product name

SAMD4B antibody

## Immunogen

[Human SAMD4B Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SAMD4B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse brain tissue, mouse kidney tissue, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of SAMD4B recombinant protein (Accession Number: XM\_011527063). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115053(SAMD4B antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using SAMD4B antibody Catalog No:115053 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).