

TXNDC9 antibody



Catalog Number: 116453

Product name

TXNDC9 antibody

Immunogen

[Human TXNDC9 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TXNDC9 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in Raji cells, HEK-293 cells, K-562 cells. Observed molecular weight by Western-blot: 25-27 kDa

Preparation

This antibody was obtained by immunization of TXNDC9 recombinant protein (Accession Number: NM_005783). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116453(TXNDC9 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:116453(TXNDC9 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)