

E2F1 antibody

Catalog Number: 110153

Product name

E2F1 antibody

Immunogen

[Human E2F1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

E2F1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse brain tissue, A431 cells, HeLa cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of E2F1 recombinant protein (Accession Number: NM_005225). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110153(E2F1 antibody) at dilution of 1:50



(under 10x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110153(E2F1 Antibody) at dilution of 1:300