

PTPRN2 Polyclonal Antibody



Catalog Number: 162253

Product name

PTPRN2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PTPRN2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 22-242 of human PTPRN2 (NP_570858.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Western blot - PTPRN2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PTPRN2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.




Immunohistochemistry - PTPRN2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using PTPRN2 antibody at dilution of 1:100 (40x lens).

PTPRN2 Polyclonal Antibody




Catalog Number: 162253




Immunohistochemistry - PTPRN2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using PTPRN2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - PTPRN2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using PTPRN2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - PTPRN2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using PTPRN2 antibody at dilution of 1:100 (40x lens).