

RSRC1 Polyclonal Antibody



Catalog Number: 165253

Product name

RSRC1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RSRC1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 117-276 of human RSRC1 (NP_001258763.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - RSRC1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RSRC1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - RSRC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using RSRC1 antibody at dilution of 1:100



(40x lens).


RSRC1 Polyclonal Antibody



Catalog Number: 165253

Immunohistochemistry - RSRC1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using RSRC1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - RSRC1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using RSRC1 antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - RSRC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using RSRC1 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - RSRC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using RSRC1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - RSRC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using RSRC1 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - RSRC1 Polyclonal Antibody

RSRC1 Polyclonal Antibody



Catalog Number: 165253

Immunohistochemistry of paraffin-embedded mouse brain using RSRC1 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - RSRC1 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RSRC1 antibody .