

TOR2A antibody



Catalog Number: 116253

Product name

TOR2A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TOR2A Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001085347). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116253(TOR2A antibody) at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using Catalog No:116253(TOR2A Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG