

# Anti-DBC-1 (3G4) Mouse antibody



Catalog Number: 168253

## Product name

Anti-DBC-1 (3G4) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

DBC-1 (3G4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human DBC1 protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal antibody in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF, IP

## Dilutions


WB: 1/500

;ICC: 1/200


## Validations



Western blot detection of DBC1 in HeLa, A549, Jurkat, C6, 3T3 and COS7 cell lysates using DBC1 mouse mAb (1:500 diluted). Predicted band size: 130KDa. Observed band size: 130KDa.



Immunocytochemistry of HeLa cells using anti-DBC1 mouse mAb diluted 1:200.



Immunoprecipitation analysis of HeLa cell lysates using DBC1 mouse mAb.