

# RHBDD1 antibody



Catalog Number: 114653

## Product name

RHBDD1 antibody

## Immunogen

[Human RHBDD1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RHBDD1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, HEK-293 cells, mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of RHBDD1 recombinant protein (Accession Number: BC027900). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114653(RHBDD1 antibody) at dilution of 1:300