

ATP5F1 antibody



Catalog Number: 108353

Product name

ATP5F1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

ATP5F1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human brain tissue, human heart tissue, human kidney tissue, human skin tissue, human testis tissue. Positive IF detected in Hela cells. Positive WB detected in mouse liver tissue, mouse heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 25kd

Validations




Preparation

This antibody was obtained by immunization of ATP5F1 recombinant protein (Accession Number: NM_001688). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108353(ATP5F1 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human liver using Catalog No:108353(ATP5F1 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF




Dilutions

ATP5F1 antibody



Catalog Number: 108353

Immunohistochemical of paraffin-embedded human liver using Catalog No:108353(ATP5F1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of Hela cells, using ATP5F1 antibody Catalog No:108353 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).