

# CNOT7 antibody



Catalog Number: 109353

## Product name

CNOT7 antibody

## Immunogen

[Human CNOT7 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CNOT7 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 33-35 kDa

## Preparation

This antibody was obtained by immunization of CNOT7 recombinant protein (Accession Number: XM\_047421731). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109353(CNOT7 antibody) at dilution of 1:500



IP Result of anti-CNOT7 (IP:Catalog No:109353, 4ug; Detection:Catalog No:109353 1:300) with mouse testis tissue lysate 4800ug.