

PCSK7 antibody



Catalog Number: 113753

Product name

PCSK7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCSK7 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of PCSK7 recombinant protein (Accession Number: NM_004716). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113753(PCSK7 antibody) at dilution of 1:50



(under 10x lens)

mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113753(PCSK7 Antibody) at dilution of 1:300