

Anti-ATG3 Rabbit antibody



Catalog Number: 166753

Product name

Anti-ATG3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATG3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ATG3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000-1/10000

IHC: 1/50


ICC/IF: 1/50

FC: 1/200

Validations



Immunocytochemistry of ATG3(green) in HeLa cells using ATG3 Rabbit pAb at dilution 1/50, and



Western blot detection of ATG3 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using ATG3 Rabbit pAb(1:1000 diluted). Predicted band size: 36kDa. Observed band size: 40kDa.