

# Anti-P44/42 MAPK(ERK1/2) Mouse antibody



Catalog Number: 169063

## Product name

Anti-P44/42 MAPK(ERK1/2) Mouse antibody

## Dilutions

1:1000-2000

## Specificity


1:50-100

Human, Rat, Mouse

## Validations

## Antibody description

P44/42 MAPK(ERK1/2) Mouse monoclonal antibody



## Preparation

Antigen: Synthetic Peptide of P44/42 MAPK (ERK1/2)

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using P44/42 MAPK (ERK1/2) Mouse mAb diluted at 1:200.

## Storage

-20°C/1 year

## Clonality


Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IHC-P



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with P44/42 MAPK(ERK1/2) Mouse mAb diluted at 1:2,000.