

WDR33 antibody



Catalog Number: 116863

Product name

WDR33 antibody

Specificity

Human; other species not tested.

Antibody description

WDR33 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 170 kDa

Preparation

This antibody was obtained by immunization of WDR33 recombinant protein (Accession Number: BC010283). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116863(WDR33 Antibody) at dilution of 1:1000



IP Result of anti-WDR33 (IP:Catalog No:116863, 4ug; Detection:Catalog No:116863 1:1000) with HeLa cells lysate 3200ug.