

# DR1 Polyclonal Antibody

Catalog Number: 161563

## Product name

DR1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to DR1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human DR1 (NP\_001929.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




## Western blot - DR1 Polyclonal Antibody

Western blot analysis of extracts of mouse spleen, using DR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



## Immunohistochemistry - DR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using DR1 antibody at dilution of 1:100 (40x lens).



# DR1 Polyclonal Antibody



Catalog Number: 161563

---

Immunofluorescence - DR1 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using DR1 antibody .